


BUDGETCARD Inc.

TT20-D: Tub Temperature Testing Card, *Stock Dolphin*

TT20-D: front



TT20-D: back

- Make sure hot water heater is set at 120 degrees.
- Fill tub with hot and cold water at the same time.
- Fill tub no higher than baby's belly button.
- Use this card and your hand when checking water temperature. If skin flushes, water is too hot.
- Put your child in the water with their back to the faucet so they cannot turn the water on.
- Never leave baby unattended.
- Empty bath as soon as you are done.

This card is a guide and is to be used in conjunction with the hand method

Please remove all guideline indicators before submitting artwork

 = **Final Cut Size: 2.125" x 3.375"**

 = **Ad Copy Area: 1.875" x .375"**
Art must be within clearances - 0.125" on all sides.